

Product Datasheet

Goat Anti-Rabbit IgG(H+L) antibody (HRP) (orb688931)

Description	Goat Anti-Rabbit IgG(H+L) antibody (HRP)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	IgG
Preservatives	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA:1:500-1:10000, WB:1:3000-1:5000, IHC:1:50-1:500, ICC:1:50-1:500
Expiration Date	12 months from date of receipt.