

Product Datasheet

HRP-conjugated Affinipure Goat Anti-Rabbit IgG(H+L) antibody (HRP) (orb688932)

Description	HRP-conjugated Affinipure Goat Anti-Rabbit IgG(H+L) antibody (HRP)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Rabbit IgG(H+L)
Preservatives	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.