



Product Datasheet

ACLY antibody (orb688949)


Description	Rabbit monoclonal antibody to ACLY
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Immunogen	A synthesized peptide derived from human ATP citrate lyase
Target	ACLY
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P53396
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HeLa Cells u...



Western blot analysis of L02 whole cell ...



Immunofluorescence analysis of HepG2 Cel...