



Product Datasheet

AURKB antibody (orb688961)


Description	Rabbit monoclonal antibody to AURKB
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Aurora B
Target	AURKB
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q96GD4
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:200, IP:1:200-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ce...



Immunoprecipitation analysis of HeLa who...



Western blot analysis of HeLa whole cell...