



Product Datasheet

PIK3CB antibody (orb688965)

Description	Rabbit monoclonal antibody to PIK3CB
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	A synthesized peptide derived from human PI3 Kinase p110 beta
Target	PIK3CB
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P42338
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of K562 who...



Western blot analysis of HeLa whole cell...