



Product Datasheet

E2F1 antibody (orb689016)


Description	Rabbit monoclonal antibody to E2F1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human E2F1
Target	E2F1
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q01094
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:200, IF:1:20-1:200, FC:1:20-1:200
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of hela cells
U...



Western blot analysis of Jurkat whole ce...



Immunohistochemical staining of human br...