


Product Datasheet

ADAM17 antibody (orb689125)

| | |
|----------------------------|---|
| Description | Rabbit monoclonal antibody to ADAM17 |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF |
| Immunogen | A synthesized peptide derived from human ADAM17 |
| Target | ADAM17 |
| Preservatives | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Monoclonal |
| Uniprot ID | P78536 |
| Dilution Range | Recommended dilution: IF:1:20-1:200 |
| Expiration Date | 12 months from date of receipt. |



Immunofluorescence analysis of HeLa Cell...