



Product Datasheet

PDCD6IP antibody (orb689145)

Description	Rabbit monoclonal antibody to PDCD6IP
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	A synthesized peptide derived from human ALIX
Target	PDCD6IP
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q8WUM4
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of Jurkat cells...



Western blot analysis of L02 whole cell ...