

## Product Datasheet

### SARS Antibody (orb689344)

<b>Description</b>	Anti-SARS-CoV-2 Spike S2 antibody(DM25), Rabbit mAb
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Target</b>	S protein
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Dilution Range</b>	ELISA 1:5000-10000
<b>Expiration Date</b>	12 months from date of receipt.