

Product Datasheet

SARS Antibody (orb689346)

| | |
|----------------------------|--|
| Description | Anti-SARS-CoV-2 RBD antibody(DM26), Rabbit mAb |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC |
| Target | S protein RBD |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Application notes | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Isotype | Rabbit IgG |
| Clonality | Recombinant |
| Uniprot ID | PODTC2 |
| Dilution Range | ELISA 1:5000-10000; Flow Cyt 1:100 |
| Expiration Date | 12 months from date of receipt. |