

Product Datasheet

SARS Antibody (orb689354)

Description	Anti-SARS-CoV-2 RBD antibody(DM35), Rabbit mAb
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	S protein RBD
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	Rabbit IgG
Clonality	Recombinant
Uniprot ID	P0DTC2
Dilution Range	ELISA 1:5000-10000; Flow Cyt 1:100
Expiration Date	12 months from date of receipt.