

Product Datasheet

SARS Antibody (orb689361)

Description	Anti-SARS-CoV-2 Spike antibody(DM42), Rabbit mAb
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Target	S protein
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	Rabbit IgG
Clonality	Recombinant
Uniprot ID	PODTC2
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.