

Product Datasheet

Clenbuterol Antibody (orb68972)

Description	Mouse monoclonal antibody to Clenbuterol
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Clenbuterol-KLH.
Preservatives	Purified antibody in PBS containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.