


Product Datasheet

KSHV ORF8 (orb68988)

Description	Rabbit polyclonal antibody to KSHV ORF8
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Western blot analysis of uninduced BCBL1...