



Product Datasheet

S100A10/P11 Antibody (orb69014)


Description	Mouse monoclonal antibody to S100A10/P11
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Purified recombinant fragment of human P11 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	11kDa
Entrez	6281
Dilution Range	WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Immunohistochemical analysis of paraffin...



Western blot analysis of MCF-7 (lane1), ...