


Product Datasheet

Influenza B virus Nucleoprotein Antibody (orb69084)

Description	Mouse monoclonal antibody to Influenza B virus Nucleoprotein
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of Influenza B virus Nucleoprotein (strain:B/Lee/40) expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
Entrez	26824002
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of full-length rec...