


Product Datasheet


GRAP2 antibody (PE) (orb691932)

Description	Mouse monoclonal antibody against GRAP2 conjugated to PE
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	GST-fusion human GRAP2/GADS protein
Target	GRAP2
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	075791
Entrez	9402
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Separation of lymphocytes stained using

...



Flow cytometry intracellular staining pa...