

## Product Datasheet

### LLY-283 (orb691941)

<b>Description</b>	LLY-283 (LLY283) is a potent, selective inhibitor of arginine methyltransferase 5 (PRMT5) with IC50 ...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{17}H_{18}N_4O_4$
<b>MW</b>	342.36
<b>CAS Number</b>	[2040291-27-6]
<b>SMILES</b>	<chem>O1[C@H]([C@H](O)C2=CC=CC=C2)[C@@H](O)[C@@H](O)[C@@H]1N1C2C(C=C1)=C(N)N=CN=2</chem>
<b>Expiration Date</b>	12 months from date of receipt.