

## Product Datasheet

### Egfr-IN-7(TQB3804 ) (orb691994)

<b>Description</b>	TQB3804 is a selective and potent EGFR kinase inhibitor, has IC50 values of 7.92 nM,0.218 nM and 0.1...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{32}H_{41}BRN_9O_2P$
<b>MW</b>	694.6048057
<b>CAS Number</b>	[2267329-76-8]
<b>SMILES</b>	<chem>BrC1=CN=C(N=C1NC1C=CC2C(C=1P(C)(C)=O)=NC=CN=2)NC1=C(C=C(C(C)=C1)N1CCC(CC1)N1CCN(C)CC1)OC</chem>
<b>Expiration Date</b>	12 months from date of receipt.