

Product Datasheet

VU 0465350 (orb692051)

Description	VU 0465350 is a potent, specific BMP receptor activin receptor-like kinase 3 (ALK3) antaognist with ...
Storage	2 years -20°C Powder, 2 weeks 4°C in DMSO, 6 months -80°C in DMSO
Note	For research use only
Formula	$C_{19}H_{18}N_4O$
Purity	> 98%
MW	318.38
CAS Number	[1414843-61-0]
SMILES	<chem>C12=NC=C(C3=CN=C3)N1C=CC(C1=CC=C(OC(C)C)C=C1)=C2</chem>
Expiration Date	12 months from date of receipt.