

## Product Datasheet

### TH-257 (orb692057)

<b>Description</b>	TH-257 (TH257) is a potent, selective, allosteric inhibitor of LIMK1 and LIMK2 with IC50 of 84 nM fo...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{24}H_{26}N_2O_3S$
<b>MW</b>	422.543
<b>CAS Number</b>	[2244678-29-1]
<b>SMILES</b>	<chem>C(N(CC1=CC=CC=C1)CCCC)(=O)C1=CC=C(S(NC2=CC=CC=C2)(=O)=O)C=C1</chem>
<b>Expiration Date</b>	12 months from date of receipt.