

## Product Datasheet

### ML417 (orb692122)

<b>Description</b>	ML417 is a Highly Selective D3 Dopamine Receptor Agonist. ML417 potently promotes D3R-mediated $\beta$ -ar...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{22}H_{25}N_3O_3$
<b>MW</b>	379.45
<b>CAS Number</b>	[1386162-69-1]
<b>SMILES</b>	<chem>O=C(C(N1)=CC2=C1C=CC=C2)N3CCN(CCOC4=CC=C(OC)C=C4)CC3</chem>
<b>Expiration Date</b>	12 months from date of receipt.