

Product Datasheet

ML401 (CID73169083) (orb692143)

| | |
|------------------------|---|
| Description | ML401 (CID73169083) is a potent, functional antagonist of EB1-2 (Epstein-Barr virus-induced gene 2) ... |
| Note | For research use only |
| Formula | $C_{20}H_{20}BrClN_2O$ |
| MW | 419.747 |
| CAS Number | [1597489-14-9] |
| SMILES | <chem>C1CN(CCN1CC2=CC=C(C=C2)Cl)C(=O)/C=C/C3=CC=C(C=C3)Br</chem> |
| Expiration Date | 12 months from date of receipt. |