

Product Datasheet

HighFidelity RED PCR Labeling Testkit (orb692243)

Description HighFidelity RED PCR Labeling Testkit is designed to produce randomly Cy5-, AF647-and ATTO643-modifi...

Storage store at -20 °C

Note For research use only

Application notes **Spectroscopic Propertie:** Cy5: λ_{exc} 649 nm, λ_{em} 670 nm, ϵ 250.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5), AF647: λ_{exc} 648 nm, λ_{em} 671 nm, ϵ 270.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5), ATTO643: λ_{exc} 643 nm, λ_{em} 665 nm, ϵ 150.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5).

Hazard Information Non-Toxic

Expiration Date 12 months from date of receipt.