

Product Datasheet

Other BK ELISA Kit (Ready to Use) (orb692743)

Description	Competitive Inhibition Other Bradykinin (BK) ELISA Kit Ready to Use: Pre-diluted detection reagents f...
Reactivity	All
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	8.65pg/mL
Expiration Date	Please enquire.