

Product Datasheet

Mouse INHbA ELISA Kit (Ready to Use) (orb692780)

Description	Sandwich Mouse Inhibin Beta A (INHbA) ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q04998
Sensitivity	6.22pg/mL
Entrez	16323
Expiration Date	Please enquire.