

Product Datasheet

Canine IL12B ELISA Kit (orb692816)

Description	Sandwich Canine Interleukin 12B (IL12B) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	6.65pg/mL
Expiration Date	Please enquire.