

Product Datasheet

Other LXA4 ELISA Kit (orb692832)

Description	Competitive Inhibition Other Lipoxin A4 (LXA4) ELISA Kit
Reactivity	All
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.132ng/mL
Expiration Date	Please enquire.