

Product Datasheet

Other LXA4 ELISA Kit (Ready to Use) (orb692833)

Description	Competitive Inhibition Other Lipoxin A4 (LXA4) ELISA KitReady to Use: Pre-diluted detection reagents...
Reactivity	All
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.132ng/mL
Expiration Date	Please enquire.