

Product Datasheet

Guinea pig IL2 ELISA Kit (orb692840)

Description	Sandwich Guinea pig Interleukin 2 (IL2) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	6.9pg/mL
Expiration Date	Please enquire.