

Product Datasheet

Human ATG3 ELISA Kit (orb692848)

Description	Sandwich Human Autophagy Related Protein 3 (ATG3) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	6.9pg/mL
Expiration Date	Please enquire.