

## Product Datasheet

### Human GnRHR ELISA Kit (orb692864)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Sandwich Human Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit                 |
| <b>Reactivity</b>        | Human  |
| <b>Range</b>             | 0.156-10ng/mL  |
| <b>Form/Appearance</b>   | Traditional  |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | serum, plasma, tissue homogenates and other biological fluids.                           |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates and other biological fluids.                           |
| <b>Sensitivity</b>       | 0.07ng/mL  |
| <b>Expiration Date</b>   | Please enquire.  |