

Product Datasheet

Human NAP2 ELISA Kit (orb692874)

Description	Sandwich Human Neutrophil Activating Protein 2 (NAP2) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q99733
Sensitivity	10.1pg/mL
Entrez	4676
Expiration Date	Please enquire.