

## Product Datasheet

### Human PANX1 ELISA Kit (Ready to Use) (orb692877)

<b>Description</b>	Sandwich Human Pannexin-1 (PANX1) ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q96RD7</a>
<b>Sensitivity</b>	0.060ng/mL
<b>Entrez</b>	<a href="#">24145</a>
<b>Expiration Date</b>	Please enquire.