

Product Datasheet

Monkey INS ELISA Kit (orb692890)

Description	Competitive Inhibition Monkey Insulin (INS) ELISA Kit
Reactivity	Monkey
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	0.12ng/mL
Expiration Date	Please enquire.