

## Product Datasheet

### Mouse NLRP3 ELISA Kit (orb692915)

<b>Description</b>	Sandwich Mouse NACHT, LRR and PYD Domains-Containing Protein 3(NLRP3) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q8R4B8</a>
<b>Sensitivity</b>	0.055ng/mL
<b>Entrez</b>	<a href="#">216799</a>
<b>Expiration Date</b>	Please enquire.