



Product Datasheet

HAS1 Antibody (orb69294)

Description	Mouse monoclonal antibody to HAS1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Purified recombinant fragment of human HAS1 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	65kDa
Entrez	3036
Dilution Range	ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U251 cell...



Western blot analysis of human HAS1 reco...