

## Product Datasheet

### MAPKAPK2 antibody (orb693482)

<b>Description</b>	Rabbit polyclonal antibody to MAPKAPK2
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Immunogen</b>	Synthetic Peptide (N-term)
<b>Target</b>	MAPKAPK2
<b>Form/Appearance</b>	Supplied in crude ascites with 0.09% (W/V) sodium azide.
<b>Storage</b>	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purified
<b>MW</b>	45568 Da
<b>Uniprot ID</b>	<a href="#">P49137</a>
<b>NCBI</b>	<a href="#">NP_004750.1</a> , <a href="#">NP_116584.2</a>
<b>Dilution Range</b>	WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100
<b>Expiration Date</b>	12 months from date of receipt.