



Product Datasheet

CD94 Antibody (orb69388)


Description	Mouse monoclonal antibody to CD94
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human CD94 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	21kDa
Entrez	3824
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
RAJI cells u...



Western blot
analysis of
human CD94
(AA:...



Line graph
illustrates
about the AG-
Ab r...