



Product Datasheet

Tubulin alpha 1A antibody (orb69674)


Description	Tubulin alpha 1A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q71U36
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell lysa...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...