



Product Datasheet

PYST2 antibody (orb69680)


Description	PYST2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16829
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell lysa...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...