



Product Datasheet

MPI antibody (orb69710)

Description	MPI antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P34949
Expiration Date	12 months from date of receipt.



Western blot analysis of A: 293T, B:A431...



Immunohistochemical analysis of formalin...