


Product Datasheet

OTUB1 antibody (orb69715)

Description	OTUB1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q96FW1
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:HeLa, B:HepG2...



Confocal immunofluorescence analysis of ...