

Product Datasheet

glutathione S-transferase alpha 2 antibody (orb69801)

Description glutathione S-transferase alpha 2 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, IHC, WB


Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note For research use only


Clonality Recombinant

Uniprot ID **P09210**


Expiration Date 12 months from date of receipt.




Immunohistochemical analysis of formalin...



Costained with Hoechst 33342



Confocal immunofluorescence analysis of ...



Western Blot analysis of A:HepG2 whole c...