



Product Datasheet

SEC13 protein isoform 1 antibody (orb69821)

Description	SEC13 protein isoform 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P55735
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:293T, B:A431,...



Immunohistochemical analysis of formalin...