



Product Datasheet

Edc3 antibody (orb69931)

Description	Edc3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q96F86
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:A431 whole ce...



Immunohistochemical analysis of formalin...