

Product Datasheet

RAB1A antibody (orb699461)

Description	Mouse monoclonal antibody against RAB1A conjugated to Biotin
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	IP, WB
Immunogen	Full length recombinant RAB1A protein
Target	RAB1A
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-613 was sufficient for detection of RAB1A in 10 ug (A549, MEF); 20 ug (Mouse Brain) lysates by ECL immunoblot analysis.
Isotype	IgG1
Clonality	Monoclonal
MW	23 kDa
Uniprot ID	P62820
NCBI	NP_004152.1
Entrez	5861
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.