

Product Datasheet

Chromogranin A Antibody / CHGA (orb699498)

Description

Chromogranin A (CGA) is an acidic protein which is



Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

A portion of amino acids 400-600 from the human protein was used as the immunogen for the recombinant Chromogranin A antibody.


IHC staining
of FFPE
human
pancreas
with...

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the recombinant Chromogranin A antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



Note

For research use only

Application notes

Optimal dilution of the recombinant Chromogranin A antibody should be determined by the researcher.

IHC staining
of FFPE
human
pancreas
with...

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID

P10645

Hazard

Information

This recombinant Chromogranin A antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.