


Product Datasheet

Factor XIIIa Antibody / F13A1 (orb699542)

Description	It recognizes a protein of 83kDa, which is identified as Factor XIIIa. It has been identified in pla...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	A recombinant fragment of human Factor XIIIa was used as the immunogen for the Factor XIIIa antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Factor XIIIa antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the Factor XIIIa antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P00488
Hazard Information	This Factor XIIIa antibody is available for research use only.
Dilution Range	ELISA: order Ab without BSA for coating
Expiration Date	12 months from date of receipt.



SDS-PAGE analysis of purified, BSA-free ...