


Product Datasheet

CELA3B Antibody / Elastase 3B (orb699545)

Description	This MAb recognizes a protein of ~30kDa, identified as CELA3B (Chymotrypsin like elastase family mem...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A portion of amino acids 82-238 from the human protein was used as the immunogen for the CELA3B antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the CELA3B antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the CELA3B antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG3, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P08861
Hazard Information	This CELA3B antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, ELISA: 2-4ug/ml; order Ab without BSA for coating
Expiration Date	12 months from date of receipt.



Western blot testing of human PANC-1 cel...



SDS-PAGE analysis of purified, BSA-free ...